

## **COVID-19 Antigen Rapid Test**

zugelassen bei BfArM

**Handelsname: COVID-19 Antigen Rapid Test Cassette**  
**Hersteller: Clongene Biotech Co., Ltd.**

### **Inhalt**

- 1. Bilder vom Hersteller**
- 2. Test Inhalt**
- 3. Gebrauchsanweisung**
- 4. Validierungsdaten**

## COVID-19 Antigen Rapid Test



[www.bauer-karlsruhe.de](http://www.bauer-karlsruhe.de)

## COVID-19 Antigen Rapid Test

### Test Inhalt





## COVID-19 Antigen Rapid Test

### Validierungsdaten

CE49PT0505

#### Clinical Sensitivity and Specificity Study Report

The “CLUNGENE<sup>®</sup> COVID-19 Antigen Rapid Test Cassette” manufactured by Hangzhou Clongene Biotech Co.,Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

The clinical performance of the “CLUNGENE<sup>®</sup> COVID-19 Antigen Rapid Test Cassette” was assessed at clinical sites evaluation with nasopharyngeal swabs obtained from individuals suspected of with COVID-19 infection.

The clinical evaluation including 32 clinical positive specimens from individuals were finally confirmed positive for SARS-CoV-2 virus infection by RT-PCR. 130 clinical negative specimens from individuals were finally confirmed negative for SARS-CoV-2 virus infection by RT-PCR.

#### 1. Method

Regarding the SARS-CoV-2 nucleocapsid antigen detection, testing was performed on 162 clinical nasopharyngeal swab specimens. 32 positive specimens and 130 negative specimens were compared to RT-PCR. We calculate the PPA, NPA with 95%CI.

#### 2. Sample and Collection

2.1 Sample used for RT-PCR and SARS-COV-2 Antigen: nasopharyngeal swab

2.2 Sample Collection: The samples we used for SARS-COV-2 Antigen detection were retrospective samples in the COVID-19 2020 outbreak period and high risk area. For positive sample, we included confirmed SARS-CoV-2 Virus infection person with swabs tested RT-PCR positive. For negative sample, we included the swabs from asymptomatic subjects whom live in high risk area when the epidemic was under control, with swabs tested RT-PCR negative. We blind - coded the samples for testing.

2.3 Sample Storage: Freshly collected specimens were processed and tested in one hour after specimen collection. Specimen stored at 2-8°C for no more than 24 hours. Store at -70 °C for a long time.

#### 3. Comparator method

RT-PCR

#### 4. Operators

They were trained of Operating procedures for CLUNGENE<sup>®</sup> COVID-19 Antigen Rapid Test Cassette user manual and clinical evaluation protocol.

#### 5. Enrollment criteria (inclusion/exclusion criteria)

##### 5.1 Inclusion criteria

- Individuals suspected of COVID-19 infection tested by SARS-CoV-2 RT-PCR test
- Confirmed infected by SARS-CoV-2 RT-PCR test
- Asymptomatic subjects underwent SARS-CoV-2 RT-PCR test

##### 5.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

1/2

## 6. Result

COVID-19 Antigen		RT-PCR		Total
		Positive	Negative	
CLUNGENE®	Positive	27	0	27
	Negative	5	130	135
Total		32	130	162

Positive Percent Agreement (PPA)= 84.38% (27/32), (95%CI: 68.25%~93.14%)

Negative Percent Agreement (NPA) =100% (130/130), (95%CI: 97.13%~100%)

## 7. Conclusion

The clinical research is a qualitative test comparison to evaluate the clinical use validity and group professional test applicability of the “COVID-19 Antigen Rapid Test Cassette” manufactured by Hangzhou Clongene Biotech Co., Ltd.

For COVID-19 Antigen, when compared to RT-PCR, a statistical comparison was made between the results yielding PPA of 84.38% (95% CI: 68.25%~93.14%), NPA of 100% (95% CI: 97.13%~100%).

For COVID-19 Antigen Rapid Test Cassette is intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens. There were still may appear false positive results and negative results. A false negative result can occur if the quantity of antigens for the SARS-CoV-2 virus present in the specimen is below the detection threshold of the assay, or the virus has undergone minor amino acid mutation(s) in the target epitope region recognized by the monoclonal antibodies utilized in the test. Improper specimen collection, improper specimen storage or repeated freezing and thawing of specimens also can lead to inaccurate results.

Hangzhou Clongene Biotech Co., Ltd.

**COVID-19 Antigen Rapid Test**  
**Cross-reactivity (Analytical Specificity) Study Report**

Author: Qingqing Chen

Final report date: 2020.06.17

Management of the study: Hangzhou Clongene Biotech Co., Ltd.  
R& D Department  
Quality Management department

Place of study: Hangzhou

Sponsor: Hangzhou Clongene Biotech Co., Ltd.

### Study Director Signature and Verification Dates

Study Director: Qingqing Chen

Company: Hangzhou Clongene Biotech Co., Ltd

Position: Director of R&D Department

Signature: 陈青青

Date: 2020.06.17

Study Handlers: Lulu Zhang

Signature: 张露

Date: 2020.06.17

Verifier: Zhiqiang Yin

Signature: 殷智强

Date: 2020.06.17

The study dates were as follows:

Protocol Effective Date: 2020-06-15

Test Starting Date: 2020-06-16

Test Duration Date: 2020-06-16~ 2020-06-16

Final Report Date: 2020-06-17



## Study Summary

We tested related pathogens, high prevalence disease agents and pathogenic flora that are reasonably likely to be encountered in the clinical specimen by COVID-19 Antigen Rapid Test. No cross-reactivity is observed within virus/bacteria at the certain concentration in the assay.

## 1. Purpose

The following study was done to confirm that the test does not react with related pathogens, high prevalence disease agents and pathogenic flora that are reasonably likely to be encountered in the clinical specimen.

## 2. Reference

The study was conducted according to Technical Specification of COVID-19 Antigen Rapid Test.

## 3. Materials

Three sequential batches of COVID-19 Antigen Rapid Test (Lot#: 2020050128, 2020050129, 2020050130)

## 4. Method

Some related pathogens, high prevalence disease agents and pathogenic flora that are reasonably likely to be encountered in the clinical specimen with certain concentration were tested by each batch of the COVID-19 Antigen Rapid Test. Every specimen were tested for three times.

## 5. Acceptable standard

The Virus or Bacteria culture with certain concentration in negative clinical matrix should not influence the negative results.

## 6. Result

Carrying out tests in accordance with standard operating procedures for use of the COVID-19 Antigen Rapid Test Cassette. The results as follows (“+”=Positive, “-”= Negative):

Virus/Bacteria	Concentration	2020050128	2020050129	2020050130
Influenza A (H1N1)	1×10 <sup>6</sup> PFU/mL	-	-	-
Influenza A (H3N2)	1×10 <sup>6</sup> PFU/mL	-	-	-
Influenza B (Yamagata)	1×10 <sup>6</sup> PFU/mL	-	-	-
Influenza B (Victoria)	1×10 <sup>6</sup> PFU/mL	-	-	-
Adenovirus	1×10 <sup>6</sup> PFU/mL	-	-	-
Human metapneumovirus	1×10 <sup>6</sup> PFU/mL	-	-	-
Parainfluenza virus	1×10 <sup>6</sup> PFU/mL	-	-	-
Respiratory syncytial virus	1×10 <sup>6</sup> PFU/mL	-	-	-
Streptococcus pyogenes	1×10 <sup>7</sup> CFU/mL	-	-	-
Candida albicans	1×10 <sup>7</sup> CFU/mL	-	-	-
Mycoplasma pneumoniae	1×10 <sup>7</sup> CFU/mL	-	-	-
Chlamydia pneumoniae	1×10 <sup>7</sup> CFU/mL	-	-	-

Legionella pneumophila	1×10 <sup>7</sup> CFU/mL	-	-	-
Human coronavirus 229E	1×10 <sup>6</sup> PFU/mL	-	-	-
Human coronavirus OC43	1×10 <sup>6</sup> PFU/mL	-	-	-
Human coronavirus NL63	1×10 <sup>6</sup> PFU/mL	-	-	-
Human coronavirus HKU1	1×10 <sup>6</sup> PFU/mL	-	-	-

## 7. Conclusion

Related pathogens, high prevalence disease agents and pathogenic flora that are reasonably likely to be encountered in the clinical specimen were tested. It showed that no cross-reactivity was found in these tests.

Hangzhou Clongene Biotech Co., Ltd.

## **COVID-19 Antigen Rapid Test**

### **Stability Study Report**

Author: Qingqing Chen

Final report date: 2020.07.13

Management of the study: Hangzhou Clongene Biotech Co., Ltd.  
R & D Department  
Quality Management Department

Place of study: Hangzhou

Sponsor: Hangzhou Clongene Biotech Co., Ltd.

### Study Director Signature and Verification Dates

Study Director: Qingqing Chen

Company: Hangzhou Clongene Biotech Co., Ltd

Position: Director of R&D Department

Signature: 陈青青

Date: 2020.07.13

Study Handlers: Lulu Zhang

Signature: 张璐

Date: 2020.07.13

Verifier: Zhiqiang Yin

Signature: 殷智强

Date: 2020.07.13

The study dates were as follows:

Protocol Effective Date: 2020-04-09

Test starting Date: 2020-04-10

Test duration date: 2020-04-10 ~ 2020-07-10

Final Report Date: 2020-07-13

## Study Summary

The COVID-19 Antigen Rapid Test, from batch COVG20200401, COVG20200402, COVG20200403 taken for accelerated stability study, in-use stability, transport stability. For accelerated stability study, kits were stored under 55°C for 91 days at accelerated aging conditions for a shelf-life test of 24 months, still meet the related test criterion. For In-use stability, after the first opening of the primary container under the condition: ①temperature: 15 °C , relative humidity: 45-65%; ② Temperature 30 °C , relative humidity 45-65%; ③ Temperature 40°C, relative humidity 45-65%; ④Temperature 20°C, relative humidity > 80%; ⑤ Temperature 30°C, relative humidity > 80%; ⑥Temperature 40°C, relative humidity > 80%. is ideal in 2 hours. The transport simulated for storing at -20°C and 40°C for 35 days still meet the related test criterion.

## 1. Purpose

The objective and purpose of this stability studies is to make sure the product's shelf-life, transport and storage conditions, and stability of the reagent in-use after the first opening of the primary container.

## 2. References

The study was conducted according to Technical Specification of COVID-19 Antigen Rapid Test.

## 3. Materials and Equipment

3.1 Reagents: three batches of the COVID-19 Antigen Rapid Test kits (Lot#: COVG20200401, COVG20200402, COVG20200403).

3.2 Micrometer, Timer

## 4. Method

In this test, the stability was evaluated under the following three conditions:

### 4.1 Accelerated stability test

Formula to decide the shelf -life (according to Arrhenius's equation):

$$AAF = Q_{10}^{(T_{AA} - T_{RT})/10}$$

AAF = Accelerated aging factor

$Q_{10}$  = An aging factor for 10°C increase or decrease in temperature (Using the Arrhenius equation with  $Q_{10}$  equal to 2 is a means of calculating an aging factor.)

$T_{AA}$  = Accelerated aging temperature (55°C),

$T_{RT}$  = Ambient temperature (25°C): Select a temperature that represents the actual product storage and use conditions (NOTE 1—This temperature is typically between 20 to 25°C. A temperature of 25°C is considered a conservative approach).

Where  $Q_{10} = 2$ ; ambient temperature = 25°C; accelerated aging temperature = 55°C;

$$AAF = 2.0^{(55-25)/10};$$

$$AAF = 2.0^{3.0} = 8;$$

$$AAT = 365 * 2 \text{ days} / 8; \text{ and}$$

$$AAT = 91 \text{ days at accelerated aging conditions for a shelf-life test of 24 months (real-time equivalent).}$$

Method: Take enough COVID-19 Antigen Rapid Test from the three sequential batches of the COVID-19 Antigen Rapid Test (Lot#: COVG20200401, COVG20200402, COVG20200403) and store them under 55°C for 91 days. The relative humidity is normal. The performance is evaluated at 7<sup>th</sup>, 14<sup>th</sup>, 21<sup>st</sup>, 28<sup>th</sup>, 35<sup>th</sup>, 42<sup>nd</sup>, 56<sup>th</sup>, 70<sup>th</sup>, 77<sup>th</sup>, 84<sup>th</sup>, 91<sup>st</sup> days. Observe and record down the result of stability research.

Test duration date: 2020-04-10 ~ 2020-07-10

The amount of each test: 19 copies.

Test item:

Test Item	Test Sample	Test Content
Physical property	Extraction Reagent	Testing each sample three times
Specificity	10 COVID-19 Antigen negative controls	Testing each sample once
Sensitivity	1 low positive control (L:50ng/mL), 1 middle positive control (M:1000ng/mL)	Testing each sample three times

### 4.2 In-use stability

Method: Three batches of the COVID-19 Antigen Rapid Test (Lot#: COVG20200401, COVG20200402, COVG20200403) were divided into six groups, after the aluminum foil bag was unsealed, the test reagents were placed in the following different temperature and humidity conditions in turn and the performance is tested at 1st, 2nd, 3rd hours. ①temperature: 15°C, relative humidity: 45-65%;②Temperature 30°C, relative humidity 45-65%;③Temperature 40°C, relative humidity 45-65%;④Temperature 20°C, relative humidity > 80%;⑤Temperature 30°C, relative humidity > 80%;⑥Temperature 40°C, relative humidity > 80%. Observe and record down the result of stability research.

Test duration date: 2020-04-14 ~ 2020-04-14

Test frequency: once an hour.

The amount of each test: 19 copies.

Test item:

Test Item	Test Sample	Test Content
Physical property	Extraction Reagent	Testing each sample three times
Specificity	10 COVID-19 Antigen negative controls	Testing each sample once
Sensitivity	1 low positive control(L:50ng/mL)、1 middle positive control (M:1000ng/mL)	Testing each sample three times

#### 4.3 Transport simulation

Method: The extreme shipping temperatures were simulated by storing kits at -20°C and 40°C. Take three batches of COVID-19 Antigen Rapid Test (Lot#: COVG20200401, COVG20200402, COVG20200403) by storing devices at -20°C and 40°C. The performance of the test kits is tested at 7<sup>th</sup>, 14<sup>th</sup>, 21<sup>st</sup>, 24<sup>th</sup>, 28<sup>th</sup>, 31<sup>st</sup> and 35<sup>th</sup> days.

Test duration date: 2020-04-14 ~ 2020-05-19

The amount of test: 19 copies.

Test item:

Test Item	Test Sample	Test Content
Physical property	Extraction Reagent	Testing each sample three times
Specificity	10 COVID-19 Antigen negative controls	Testing each sample once
Sensitivity	1 low positive control(L:50ng/mL)、1 middle positive control (M:1000ng/mL)	Testing each sample three times

## 5. Evaluation Criteria

- 5.1 Physical property: The surface of reagent is flat and clean, no damage, all the materials are assembled well, and no missing component. The width of the membrane shall be wider than 2.5 mm. If the liquid migration rate is not be lower than 10mm/min, then the test for qualified.
- 5.2 Specificity: All of the negative controls should be tested negative, the test is qualified, and otherwise it is unqualified.
- 5.3 Sensitivity: All the positive controls should be tested positive, the test is qualified, and otherwise it is unqualified.

## 6. Result

### 6.1 Accelerated stability study

Lot#: COVG20200401



Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
42 <sup>nd</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
56 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
70 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
77 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
84 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
91 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
42 <sup>nd</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
56 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
70 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
77 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
84 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
91 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
42 <sup>nd</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
56 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
70 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

77 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
84 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
91 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

## 6.2 In-use stability

### 6.2.1 Store kits after open the pouch at temperature: 15°C, relative humidity: 45-65%

Lot#: COVG20200401

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

### 6.2.2 Store kits after open the pouch at temperature: 30°C, relative humidity: 45-65%

Lot#: COVG20200401

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

### 6.2.3 Store kits after open the pouch at temperature: 40°C, relative humidity: 45-65%

Lot#: COVG20200401

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(9/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(9/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

### 6.2.4 Store kits after open the pouch at temperature: 20°C, relative humidity > 80%

Lot#: COVG20200401

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
----------------------	-----------	----------	--------	--------

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

### 6.2.5 Store kits after open the pouch at temperature: 30°C, relative humidity > 80%

Lot#: COVG20200401

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

### 6.2.6 Store kits after open the pouch at temperature: 40°C, relative humidity > 80%

Lot#: COVG20200401

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(9/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
1 <sup>st</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
2 <sup>nd</sup> hour	Qualified	-(10/10)	+(3/3)	+(3/3)
3 <sup>rd</sup> hour	Qualified	-(9/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Notes: The positive results of the test cassettes are not obvious relatively (The shade of color in the test region (T) becomes shallow.) when have taken for 3 hours in the abnormal conditions.

### 6.3 Transport simulation

#### Store at -20°C for transport simulation

Lot#: COVG20200401

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
24 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
31 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
24 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
31 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: "+"=Positive "-"= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
24 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
31 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: "+"=Positive "-"= Negative

**Store at 40°C for transport simulation**

Lot#: COVG20200401

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
24 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
31 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: "+"=Positive "-"= Negative

Lot#: COVG20200402

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

24 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
31 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

Lot#: COVG20200403

Time \ Results	Physical property	Specificity	Sensitivity	
			L:50ng/mL	M:1000ng/mL
7 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
14 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
21 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
24 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
28 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
31 <sup>st</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)
35 <sup>th</sup> day	Qualified	-(10/10)	+(3/3)	+(3/3)

Result specification: “+”=Positive “-”= Negative

## 7. Conclusion

The COVID-19 Antigen Rapid Test, from batch COVG20200401, COVG20200402, COVG20200403 stored under 55°C for 91 days at accelerated aging conditions for a shelf-life test of 24 months. They still meet the related test criterion. It should be reliable to define the the expiry date of the reagents under the temperature range is 4°C~30°C for 24 months.

The stability of the reagent from batch COVG20200401, COVG20200402, COVG20200403 in-use after the first opening of the primary container for 2 hours is ideal in condition: ①temperature: 15°C, relative humidity: 45-65%;②Temperature 30°C, relative humidity 45-65%;③Temperature 40°C, relative humidity 45-65%;④Temperature 20°C, relative humidity > 80%;⑤Temperature 30°C, relative humidity > 80%;⑥ Temperature 40°C, relative humidity > 80%. For better performance, we suggested tested within 1 hour after removal from pouch specially if the room temperature is more than 30°C and in high humidity environment. The transport simulated for storing at -20°C and 40°C for 35 days still meet the related test criterion. It should be reliable to say the transport stability is good.